



SIRCAL

TECHNICAL BULLETIN 11.1

How do you make sure you are ordering the right purifier tubes for your Sircal Rare Gas Purifier?

All you need to do is ...

1. Check the serial number (S/N) of your Sircal purifier.
2. A 6-digit S/N system would be the original purifier which requires:
 - i. Titanium tube P4013
 - ii. Copper oxide tube P4014
 - iii. Molecular sieve tube P4015
3. A 10-digit S/N system would be the MP-2000 purifier which requires:
 - i. Titanium tube P4003
 - ii. Copper oxide tube P4004
 - iii. Molecular sieve tube P4005

The life span of the purifier tubes can be calculated by the method explained in the instruction manual.

“As a rule of thumb the Titanium and Copper Oxide Tubes should be replaced every 2 -3 years of use.”


The replacement procedure is very simple and you can manage it easily. Please refer to the instruction manual for a detailed procedure for replacement of your purifier tubes.

For further information, please contact us now:

Sircal Instruments (UK) Limited

Tel: +44 (0) 1342 335 309
Fax: +44 (0) 1342 323 600
Email: sales@sircal.co.uk
www.sircal.co.uk

Purifier Tubes for Original Purifier (with 6-digit serial nos.)



Purifier Tubes for MP-2000 (with 10-digit serial nos.)

